

Tunable Operation

Many lasers can operate at more than one wavelength. Argon and krypton lasers can operate at discrete wavelengths ranging from the ultraviolet to the near infrared. Dye lasers can be continuously tuned over a spectrum of wavelengths determined by the fluorescence bandwidths of the specific dyes (typically about 150 nm). Alexandrite and titanium sapphire lasers can be tuned continuously over specific spectral regions.

To create a tunable laser, the cavity coatings must be sufficiently broadband to accommodate the entire tuning range, and a variable-wavelength tuning element must be introduced into the cavity, either between the cavity optics or replacing the high-reflecting optic, to introduce loss at undesired wavelengths.

Three tuning mechanisms are in general use: Littrow prisms, diffraction gratings, and birefringent filters. Littrow prisms (see figure 10.16) and their close relative, the full-dispersing prism, are used extensively with gas lasers that operate at discrete wavelengths. In its simplest form, the Littrow prism is a 30-60-90-degree prism with the surface opposite the 60-degree angle coated with a broadband high-reflecting coating. The prism is oriented so that the desired wavelength is reflected back along the optical axis, and the other wavelengths are dispersed off axis. By rotating the prism the retroreflected wavelength can be changed. In laser applications, the prism replaces the high-reflecting mirror, and the prism's angles are altered (typically to 34, 56, and 90 degrees) to minimize intracavity losses by having the beam enter the prism exactly at Brewster's angle. For higher-power lasers which require greater dispersion to separate closely spaced lines, the Littrow prism can be replaced by a full-dispersing prism coupled with a high reflecting mirror.

Gratings are used for laser systems that require a higher degree of dispersion than that of a full-dispersing prism.

Birefringent filters have come into general use for continuously tunable dye and Ti:Sapphire lasers, since they introduce significantly lower loss than do gratings. The filter is made from a thin, crystalline-quartz plate with its fast axis oriented in the plane of the plate. The filter, placed at Brewster's angle in the laser beam, acts like a weak etalon with a free spectral range wider than the gain curve of the lasing medium. Rotating the filter around the normal to its face shifts the transmission bands, tuning the laser. Since there are no coatings and the filter is at Brewster's angle (thereby polarizing the laser), there are no inherent cavity reflection losses at the peak of the transmission band. A single filter does not have as significant a line-narrowing effect as does a grating, but this can be overcome by stacking multiple filter plates together, with each successive plate having a smaller free spectral range.




Figure 10.16 Littrow prism used to select a single wavelength



Prism-tunable-ion laser